3M Self-Ligating Bracket Overview





3M Self Ligating Bracket Overview

- History of Self-Ligation
- Benefits of Self-Ligation
 - Hygiene
 - Chair Time Savings
 - Frictional Resistance
 - Biological Resistance
 - Lighter Forces
- What is Self-Ligation?
- Self-Ligating Designs
 - Past Challenges with Self-Ligation
 - Bracket Designs
 - 3M Self-Ligation



History of Self-Ligation





History of Self-Ligation

- The self-ligation concept is <u>not</u> new
 - The first self-ligating concepts were developed in 1935 and the brackets have continued to evolve over time
- Today, self-ligating brackets represent the fastest growing orthodontic bracket segment



History of Self-Ligation

1935	Russell Lock
1971	Wildman (Edgelok); Button capped bracket opened vertically first passive slot
1975	Hanson: Speed™ (Strite); Narrow single design, difficult to control rotations and finish, technique sensitive
1980	Forestadent Mobil-Lock; Common mechanism failure
1986	Pletcher (Activa) A' Co; Difficult mechanism to manage
1995	Time™ (Adenta/American Ortho); Single wing design, non-mesh base, large bracket
1995	TwinLock (Wildman) Ormco; Precursor to the eventual Damon™ III design
1996,1999	Damon™ (A Co.); Single wing design, large external sliding door, multiple malfunctions
2000	Damon™ II (Ormco); Single wing design, reduced size from original same issues
2000	In-Ovation® (GAC); Single bracket, large size frequent clip failure
2003/2004	In-Ovation-R™ (GAC) 2003/2004; Single more functional clip closed better, reduced size
2004	SmartClip™ Brackets (3M Unitek); First true twin-wing design, novel ligating mechanism, high clip forces
2004	Damon III™ (Ormco); 2004 Single bracket door shut better, rhomboid angulation, improved mechanism
2004	Time 2 (Adenta/American); Single bracket, clip moves easier
2006	Damon™ Mx (Ormco); Combination metal/polycarbonate material for aesthetics
2006	SmartClip™ Brackets (3M Unitek); Improved clip forces
2007	Clarity [™] SL Brackets (3M Unitek); Ceramic version of SmartClip Bracket, further improved clip forces
2009	SmartClip™ SL3 Brackets (3M Unitek); Further clip force reduction

All trademarks are property of their respective holders.




Benefits of Self-Ligation





Improved Hygiene

- Ligatures can retain food particles around the bracket
 - This makes it more difficult for patients to keep teeth clean and maintain effective oral hygiene
- Ligatures are made of an elastomeric material that stains over time from food products
- Self-ligating brackets reduce the need to use ligatures, thus providing improved oral hygiene







Chair-Time Savings

Doctors save time with self-ligation

- Archwire changes are shorter because they do not have to remove and re-tie ligatures (See wire comparison videos from Dr. Sondhi on 3MUnitek.com)
- 2. Chair-time may be reduced as patients no longer debate about what color of AlastiK™ Ligature to choose





What is Friction?

Friction

- The resistance to <u>motion</u> when one object moves against another
- Examples include the archwire sliding in the bracket slot, and, brackets sliding on the archwire




Frictional Resistance

Ligated Appliances

- Must secure/ligate the archwire
- This imparts a frictional resistance with AlastiK[™] and/or SS Ligatures



- Alternative for archwire ligation
- Reduces or eliminates frictional resistance
 - Active Mechanism
 - Passive Mechanism








Frictional Resistance

- Ligated brackets have more frictional resistance than self-ligating brackets due to the addition of ligatures
- Imagine driving your car with the emergency brake on ...
 - The ligatures act as the brake
 - Driving with the emergency brake engaged requires more force to move the car
- With ligated systems, doctors need to use higher forces to overcome the frictional resistance generated by the ligature






Biological Resistance


Self-Ligating appliances move teeth with light forces

- Teeth move within the bone more efficiently when lighter forces are employed
- Heavier forces can cause the periodontal ligament (cells) to react in such a way that it restricts tooth movement





Lighter Forces





Benefits of Light Forces

- Light forces optimize cell stimulation which improves tooth movement (biomechanics)
- Light forces are achieved through the interaction of the brackets' Rx, design (twin bracket with clips) and the archwire
- The lighter forces and reduced friction offered by self-ligation may allow teeth to level and align faster
 - This may result in faster overall treatment times



SmartClip[™] SL3 Self-Ligating Appliances Friction in Treatment

- Light force initial archwires ensure patient comfort and assist in improved physiological response to accelerate leveling and aligning
- "Active on Demand" feature gives the ability to apply friction where needed with the use of ligature ties
 - Aids in effective anchorage, detailing and finishing





Biomechanics and Light Forces in relation to the MBT[™] Versatile+ Appliance System

- The MBT[™] Appliance System's reduced angulations result in lighter forces to achieve the same mechanics
- The clips provide lower friction which supports the light force mechanics of the MBT Appliance System
- Lighter forces improve cell stimulation which creates more efficient tooth movement and sliding mechanics





Biomechanics and Light Forces

- Compatible with multiple treatment disciplines without major treatment changes
 - Roth*
 - Ricketts*
- Integration of Variable Prescription Orthodontics
 - Patient Centered Treatment Solution



^{*} No endorsement by the Doctor is implied

Self-Ligating Bracket Designs





Self-Ligating Bracket Designs

- Two types of self-ligating brackets
 - Passive
 - Active
- Passive
 - No frictional force on the archwire, regardless of archwire size or material
 - Archwire remains passive within the bracket slot regardless of its size
- Active
 - With a large archwire, ligating mechanism imparts a frictional force on the archwire



Self-Ligating Bracket Designs















SmartClip™ SL3 Bracket -- Key Features

- No moving parts
- "Active on demand"
- Pre-programmed release for excessive force
- Most like Edgewise bracket in appearance
- Tooth-specific bracket bases
- Easy AlastiK[™] Ligature placement
- APC[™] PLUS and APC[™] II
 Adhesive Coated System
- Treatment Philosophy of the MBT[™] Versatile+ Appliance System






Combining the Best of Clarity™ & SmartClip™ Brackets



SELF-LIGATING APPLIANCE SYSTEM





SmartClip[™] SL3 Bracket

- Unique self-ligating technology
- Fast growing self-ligating bracket

Aesthetics and Self-ligation, without Compromise.



Ceramic technologies