

# MonoCem<sup>®</sup>


SELF-ADHESIVE RESIN CEMENT




Twist to  
remove cap



Push here to bleed tube<sup>†</sup>  
before putting on mixing tip.  
Express 1-2mm onto pad  
<sup>†</sup>First use only




**1. Prepare Tooth**  
rinse and leave moist



**2. Restoration Preparation**  
place MonoCem directly into crown  
discard mixing tip and recap syringe




**3. Restoration Seating**  
place crown onto preparation  
FLASH cure excess 1-2 seconds



**4. Remove Excess Cement**  
use a suitable instrument and  
gently remove excess



**5. PFM Non-Disturb Phase**  
maintain positive pressure for 2-1/2 minutes  
then light cure margins for 20-40\* seconds  
with complete auto cure in 3 minutes  
REMOVE any additional cement



**6. Clear and Clean Contacts**  
floss up/down and pull  
through to clean margins

\*Light cure times vary. All ceramic, crown and inlay/onlay  
cementation light cure for 20-30 seconds. Post and core  
cementation light cure for 40 seconds. Refer to MONOCERM  
DIRECTIONS FOR USE for more detailed instructions.

CEMENTS



# Complementary Products

## OneGloss®



OneGloss is a one-step finisher and polisher that provides the ideal finish for all types of composite restorations simply by alternating the contact pressure as you go from finishing to polishing.

- Aluminum oxide polishers are disposable
- Available in three different shapes
- Durable enough to withstand polishing an entire arch
- Ideal for cement removal



## Super-Snap®



Super-Snap is designed for easy, fast and safe contouring, finishing and polishing of all microfilled and hybrid composites.

- Elastic silicone shank mount for easy placement of the disk onto the mandrel
- Manufactured without a metal center
- Very flexible, ultra-thin disks allow easy access to interproximal areas
- Two sizes – four grits (silicon carbide and aluminum oxide)



ABRASIVES



For more information on these and many other great products visit [www.shofu.com](http://www.shofu.com) or call **800.827.4637**

SN77-0811

Shofu Dental Corporation • San Marcos, CA