



# SNOOP™

CARIES DETECTING DYE



Dark blue color provides strong contrast with dentin and the pulp

Allows accurate detection of infected dentin

Increased reliability, no error

Cariou dentin is made up of two distinct layers:

- The layer of outer infected dentin is soft, discolored, non-vital, nonsensitive, cannot remineralize, should be removed and is stained by Snoop in 10 seconds.
- The layer of inner uninfected (affected) dentin is not infused with bacteria, is capable of remineralizing, should not be removed and is not stained by Snoop in 10 seconds.

Snoop™ distinguishes between outer infected and inner affected dentin in 10 seconds. Snoop identifies infected carious dentin and helps the practitioner preserve vital dentin that should not be removed. The caries detecting dye stains the denatured collagen that is only present in the outer infected dentin. This is an important tool for conservative dentistry.

Photos courtesy of Dr. C.H. Pameijer



Shows obvious occlusal caries.



After removing obvious caries, apply SNOOP for 10 seconds.



Rinse and remove only the stained infected dentin.



Re-apply SNOOP and rinse. No further staining means no remaining infected dentin.



SNOOP

12mL bottle

# DENTIN DESENSITIZER

FOR APPLICATION TO ALL DENTIN SURFACES

Dentin Desensitizer contains 5% glutaraldehyde in water with fluoride added to enhance stability. It is compatible with adhesives and composites as well as traditional cements, and it does not interfere with bonding.



DES

12mL bottle