

CYCLIC LOADING (VIBRATION) ACCELERATES TOOTH MOVEMENT IN ORTHODONTIC PATIENTS: A DOUBLE-BLIND, RANDOMIZED CONTROLLED TRIAL

Dubravko Pavlin DMD, MSD, Ph.D., Ravikumar Anthony MSD, Vishnu Raj DDS, MS, Peter T. Gakunga DDS, Ph.D.
Semin Orthod, <http://dx.doi.org/10.1053/j.sodo.2015.06.005>

PURPOSE

- To assess the effect of a defined low-level cyclic loading on the rate of orthodontic tooth movement.

METHODS

- Parallel, double-blind, prospective, randomized, controlled trial.
- Enrolled 45 orthodontic patients (age range 12-40 yrs) with fixed appliances and randomized into two groups.
 - AcceleDent[®] Group (N=22)
 - Sham Control Group (N=23)
- Patients underwent extraction of maxillary first premolars with maximum maxillary posterior anchorage and at least 3mm of extraction space after initial alignment.
- Cyclic loading was applied to the vibration group for 20min/day using the AcceleDent device, which delivered a force of 0.25N (25g) at a frequency of 30Hz.
- Control group was assigned to the same protocol, but the device could not be activated to vibrate.
- Separate canine retraction was performed on a 0.018 in stainless steel archwire and enmasse retraction with a 0.019 x 0.025 SS archwire.
- Average monthly rate of maxillary canine retraction into an extraction space was analyzed.

RESULTS


- The mean rate of movement was significantly higher for the AcceleDent Group versus the Sham Control Group (Figure 1).
 - 1.16 mm/month vs 0.79mm/month (p=0.05)
- Study outcomes indicated that AcceleDent is safe and convenient for patients' daily use.

AUTHOR CONCLUSIONS

- These results showed that low-level cyclic loading of 0.25N at 30Hz increases the rate of tooth movement when applied as an adjunct to orthodontic treatment.

ACCELERATED TOOTH MOVEMENT

CYCLIC LOADING (VIBRATION) ACCELERATES TOOTH MOVEMENT IN ORTHODONTIC PATIENTS:
A DOUBLE-BLIND, RANDOMIZED CONTROLLED TRIAL



✓ KEY POINT

This peer-reviewed randomized controlled trial demonstrates that AcceleDent accelerates the rate of tooth movement during orthodontic treatment.

acceleddent.com

OrthoAccel Technologies, Inc.
1-866-866-4919
sales@orthoaccel.com

OA-L49025. 07/15

