

## EFFECT OF MECHANICAL VIBRATION ON RESISTANCE TO SLIDING IN THE FIXED ORTHODONTIC APPLIANCE


*Kennedy D, Liu D. Electronic poster session presented at: 2014 Annual Session of the American Association of Orthodontists; 2014 Apr 25-29; New Orleans, LA.*

### PURPOSE

- To test the effect of vibration on the sliding resistance (i.e. friction and binding) of a fixed orthodontic appliance system.

### METHODS

- Engineering bench test setup that simulated tooth structure with PDL layer
- Measured static and dynamic friction between 0.022" x 0.028" SS bracket and 0.016" x 0.025" NiTi wire with and without AcceleDent
  - Pulled wire at 5 mm/min for 7 mm



### RESULTS


- With AcceleDent, static and dynamic friction between orthodontic wire and brackets in a fixed appliance system were reduced by 8.5% and 22.3% ( $p < 0.05$ ), respectively.
- Range of dynamic friction at 5 mm:
  - Control: 267-478g
  - AcceleDent: 210-345g

**TREATMENT EFFICIENCY****EFFECT OF MECHANICAL VIBRATION ON RESISTANCE TO SLIDING IN THE FIXED ORTHODONTIC APPLIANCE***(Continued)***AUTHOR CONCLUSIONS**

- Our data show that vibration (AcceleDent) reduces static and dynamic frictions between orthodontic wire and brackets in fixed appliance system by about 8.5% and 22.26%, respectively.

## ACCELEDENT® AND EFFICIENCY OF FIXED APPLIANCES

AcceleDent may enhance efficiency of fixed mechanics by reducing friction



With AcceleDent, static and dynamic friction between orthodontic wire and brackets in a fixed appliance system were reduced by 8.5% and 22.3% ( $p < 0.05$ ), respectively.\*

\*Liu D, Kennedy D. Effect of Mechanical Vibration on Resistance to Sliding in the Fixed Orthodontic Appliance. Electronic poster session presented at: 2014 Annual Session of the American Association of Orthodontists; 2014 Apr 25-29; New Orleans, LA.

### ✓ KEY POINT

AcceleDent may enhance the efficiency of fixed mechanics by reducing friction.

[acceledent.com](http://acceledent.com)

OrthoAccel Technologies, Inc.  
1.866.866.4919  
sales@orthoaccel.com

