



Sparkle® Disposable Air Water Syringe Tips

**REDUCE PATIENT
RISK**
**INCREASE WORKFLOW
PRODUCTIVITY**

- ✓ **QUAD CORE® Air Channel** ensures air/water separation.
- ✓ **NO ADAPTER REQUIRED.**
- ✓ Beveled edge provides increased patient comfort.
- ✓ Single-use eliminates cleaning and cross-contamination risks.



WHY CHOOSE SPARKLE SINGLE-USE TIPS?

Sparkle Air/Water tips are: simple to connect and lock securely; ensure distinct air/water separation for every procedure; constructed with patient comfort and instrument longevity in mind; and are disposable. Traditional metal air/water syringe tips **cannot be reliably sterilized** for two main reasons: first, the inside channel is hard to reach and cannot be visually inspected. Second, researchers have found that the irregular inner surfaces of metal syringe tips allows for the accumulation of deposits and possibly bioburden.¹ If clinicians cannot remove all foreign matter from an instrument, **it is not a viable candidate for sterilization and reuse.** Eliminate the risk associated with reusable syringe tips with Sparkle Disposable Air/Water Syringe Tips.


- Tip connects directly to your system with no adapter required.
- Rigid construction will not bend during cheek retraction.
- Translucent material allows visualization of water flow and/or blockages.
- Disposable tips save you time as compared to metal syringe tips which require cleaning, disinfection and sterilization.

Available in 6 colors

**ROUNDED
NON-SCRATCH END**
Prevents O-ring Damage



LOCKING INDICATOR & GROOVE
Ensures Tip is Securely Installed



QUAD CORE® Air Channel
Ensures Air/Water Separation



BEVELED EDGE
Provides Increased
Patient Comfort



Ordering Information

Description	Quantity	Ref. #
Clear	250/Bag	BCSAWSCL
White	250/Bag	BCSAWSWH
Assorted	250/Bag	BCSAWS
Assorted	1500/Bag	BCSAWSB

Assorted Includes: Blue, Purple, Yellow, Amber

ⓧ Ⓢ = NOT MADE WITH NATURAL RUBBER LATEX



1. Puttaiah, P., Cottone, J., Guildersleeve, J., Azmoudeh, A., Tenney, J. "Rational for Using Single-Use Disposable Air/Water Syringe Tips." Compendium, November 1999, Vol. 20 No.11:1-7.