

# Duo-Check® Sterilization Pouches

Description	Quantity	Ref. #
2" x 8" / 5 x 20 cm	200/Box	SCZ
2.25" x 4" / 6 x 10 cm	200/Box	SCXX
2.75" x 9" / 7 x 23 cm	200/Box	SCX
3.25" x 9" / 8 x 23 cm	200/Box	SCQ
3.5" x 5.25" / 9 x 13 cm	200/Box	SCXS
3.5" x 9" / 9 x 23 cm	200/Box	SCS
3.5" x 9" / 9 x 23 cm	500/Box	SCS5
4.25" x 11" / 11 x 28 cm	200/Box	SCB1X
5.25" x 10" / 13 x 25 cm	200/Box	SCM
7.5" x 13" / 19 x 33 cm	200/Box	SCL
10" x 15" / 25 x 38 cm	100/Box	SCL1015
12" x 15" / 30 x 38 cm	100/Box	SCL1215
12" x 18" / 30 x 46 cm	100/Box	SCL1218

## Pink with a Purpose™ External Process Indicator

Description	Quantity	Ref. #
3.5" x 9" / 9 x 23 cm	200/Box	SCSPP
5.25" x 10" / 13 x 25 cm	200/Box	SCMPP



**1 Secure Seal**  
Duo Check offers extra wide chevron protection against instrument breaching, while Sure-Check offers advanced triple rail seal protection.

**2 Medical Grade Paper**  
Wet strength treated, virgin medical grade paper that meets or exceeds medical pouch standards.

**3 Blue Tinted Transparent Film**  
Facilitates both easy instrument viewing and identification of punctures or tears in film.

**4 Wide Self-Seal Adhesive Strip**  
Easily and securely seals pouch, providing reliable seal until pouch is opened.

**5 Perforated Fold**  
Facilitates easy, accurate folding of adhesive strip to form airtight seal until pouch is opened.

**6 Complete Paper-Film Separation**  
Consistently clean separation



(peeling) of paper and film is due to high grade materials and adhesives.

**7 Thumb Notch**  
Assures easy opening of pouch.

**8 Six Corner Tack Seals**  
Prevents corners from curling after sterilization. Tacks resist collection of dust and contaminants near seal.

**9 Process Indicators**  
Internal and external indicators change color with sterilization heat, are easily visible, and are compatible with steam (indicators turn light brown) and EtO (indicators turn dark brown) sterilization agents.

