

Rapicide® OPA/28

HIGH-LEVEL DISINFECTANT



When safety counts and minutes matter.

Ensuring a safe and healthy environment for patients and staff just got easier. Rapicide OPA/28 provides a superior way to reprocess heat-sensitive, semi-critical medical devices that are unsuitable for sterilization.

NEW! Less Odor Detection*



NEW
SANI HLD HOLDER:
• Easy - Drain Inner Tray
• Handy Built-in Calendar Reminds You When to Change Solution.

*Compared to glutaraldehyde



crosstex.com | 888.276.7783



PROTECTS®

The Rapicide OPA/28 Advantage

Ready to Use and Lasts Longer

It's important for healthcare providers to manage budgets carefully while still delivering high-quality patient care. Rapicide OPA/28 has the fastest manual High Level Disinfection time of any product in the market. This creates faster turn-around times and provides potential cost saving opportunities by creating better efficiencies.

Rapicide OPA/28:

- Use for up to 28 days
- 2-year unopened shelf life
- 75 days once opened



Highly Effective

- Effectively Inactivates: TB, HIV2, MRSA, VRE, Hepatitis Viruses
- Non-corrosive to metals

When Safety Counts

- Has less odor than glutaraldehyde*
- Because of the low percentage of OPA and low vapor pressure irritation to the eyes, nose and throat are significantly reduced when compared to glutaraldehyde**

*CDC – Guidelines for Disinfection and Sterilization in Healthcare Facilities, 2008 – reference www.cdc.gov for additional information

**EPA – Reducing Ethylene Oxide and Glutaraldehyde Use – reference www.epa.gov

† Minimum recommended concentration

More Efficient

- 10-minute manual high-level manual disinfection time...faster than any OPA or glutaraldehyde on the market today
- The faster manual disinfection time provides one extra disinfection cycle for every five instruments reprocessed allowing you to increase operational efficiencies
- No activation required

BY THE NUMBERS

10 MINUTE manual high-level disinfection time at 20°C minimum

28 DAY reuse period



Ordering Information

Description	Quantity	Ref.#
Rapicide OPA/28 High-Level Disinfectant	1-gallon bottle	ML02-0127
Rapicide OPA/28 Test Strips	50 strips/bottle	ML02-0137
Sani HLD Holder	Holds 1 gallon of liquid	JHOLD

Test Strips Provide MRC⁺ Verification

Rapicide OPA/28 Test Strips are only compatible with our disinfectant, and make measuring the concentration:

- User friendly - strips feature a large, easy-to-read indicator pad
- Intuitive - color chart on bottle provides a pass/fail reference
- Fast - 3-second dip time and 90-second test strip read time

