Biological Monitoring Service by Crosstex® BECAUSE YOU CAN'T SEE STERILE.




Sound sterility practices are vital to protecting both patients and clinicians. According to the United States Centers for Disease Control and Prevention Guidelines for Infection Control 2003, chemical and biological monitoring of all sterilizers is a necessary part of effective infection control protocol. To help dental practices meet their monitoring needs, Crosstex International offers a complete sterility assurance program.

EVERY LOAD. Sterilization Pouches

Sure-Check® Sterilization Pouches

In-office chemical indicators feature internal/external multi-parameter indicators that change color only when all three criteria for sterilization have been met... Time/Temperature/Steam.

• For use in steam and EtO sterilizers.



Patent Pending

Sure-Check®Strips

Multi-parameter internal chemical indicators for use in all pouches (except for Sure-Check and dry heat pouches) and cassettes, indicates all conditions for sterilization have been met.

• For use in steam and chemical vapor sterilizers.



UNPROCESSED





EVERY DAY. Integrators

Steam*Plus*™ Class 5 Integrators

In-office chemical indicators for use with steam sterilization cycles to provide a distinct pass/fail result, allowing for the release of all non-implant loads before the results of the spore test are received.

• Immediate readout integrator for use in gravity, prevacuum, and flash cycles.

EVERY WEEK. Biological Monitoring System

ConFirm® 10 In-Office or Mail-In Sterilizer Monitoring Service

Meet CDC guidelines for weekly biological monitoring in less than half the time...final results in just 10 hours when using in-office ConFirm® 10 Incubator and Indicators.

Or, use the most comprehensive mail-in service in healthcare with results available in just 24 hours!

- Compatible with dry heat, chemical, EtO and steam sterilizers.
- Test results in as little as 24 hours for steam, 72 hours for chemical vapor, 7 days for dry heat and EtO upon receipt of spore strips in our laboratory.







ORDER FORM

Or	NUER FURN	/	
 Sure-Check® Sterilization Pouce For use in steam sterilizers and EtO sterilizers. When used in conjunction with weekly biological maprovide a high level of sterility assurance. US Food and Drug Administration cleared. 		QTY. 200/box QTY. 100/box SCXX2 2.25" x 4" SCL10152 10" x 1 SCX2 2.75" x 9" SCL12152 12" x 1 SCXS2 3.5" x 5.25" SCL12182 12" x 1 SCS2 3.5" x 9" QTY. SCW2 5.25" x 6.5" QTY. SCM2 5.25" x 10" SCS25 3.5" x 9" SCL2 7.5" x 13" SCL5152 5.25" x 15" SCL8162 8" x 16"	15"
Sure-Check® Strips • For use in steam and chemical vapor sterilizers. • US Food and Drug Administration cleared.		QTY. 100/bag SCK	11
SteamPlus™ Class 5 Integrators Immediate readout integrator for use in gravity, process not replace the routine use of a biological indi US Food and Drug Administration cleared.	evacuum, and flash cycles.	STEAMPlus STERILIZATION INTEGRATOR All parameters and when dark like prices for the last ANT CONTROL STEAM CASS 5650 110-1200 2 (LIMIT 350 TW CONTROL STEAM) STEAMPlus STERILIZATION INTEGRATOR All parameters are when dark like prices for bine ANT CONTROL STEAM CASS 5651 100 STEAMPlus STERILIZATION INTEGRATOR STEAMPlus STERILIZATION INTEGRATOR STEAMPlus STERILIZATION INTEGRATOR All parameters are when dark like prices for bine ANT CONTROL STEAM CASS 5651 100 STEAMPlus STERILIZATION INTEGRATOR All parameters are when dark like prices for bine ANT CONTROL STEAM CASS 5651 100 STEAMPlus STERILIZATION INTEGRATOR All parameters are when dark like prices for bine ANT CONTROL STEAM CASS 5651 100 STEAMPlus STERILIZATION INTEGRATOR All parameters are when dark like prices for bine ANT CONTROL STEAM CASS 5651 100 STEAMPlus STERILIZATION INTEGRATOR AND CONTROL STEAMPLE STERILIZATION INTEGRATOR STEAMPLE STEAMPLE STERIL	
 Compatible with dry heat, chemical, EtO and steam Your choice of immediate test failure notification by follow-up phone call within 24 hours from one of ou All test results are faxed or emailed quarterly and and Check results online 24 hours a day on our secure, or complimentary reminder service provides either an notification of a missed test date. Third party documentation of test results provides proof of testing must be provided. One of the most cost-effective methods of biological Highly resistant dual-species bacterial spores provided festing for an effective sterilization cycle. US Food and Drug Administration cleared. 	sterilizers. If ax or email with a r customer care consultants. Inually. Confidential website. If ax or email alert a tangible record when monitoring. The the most reliable method	QTY. ConFirm® Premium Mail-in (Postage Included CST060 6 Tests 3 Strip Test CST120 12 Tests 3 Strip Test CST480 48 Tests 3 Strip Test ConFirm® Value Mail-in (Postage Not Included CVT120 12 Tests 2 Strip Test CVT520 52 Tests 2 Strip Test Test strips cannot be split between two sterilizers.	
ConFirm® 10 In-Office Biologica Final results in just 10 hours! • Easy-to-use and fast! • US Food and Drug Administration cleared.	al Monitoring Sys	ConFirm 10° In-Office System: C10SK ConFirm 10 Digital Incubator (Record Keeper, Manual, Specify Voltage) C10B125 25/Box ConFirm 10 Biological Indic	cators
Customer Name:		Email:	
Office realite:	FIIONE:	EIIIaII;	

Fax Orders to: 303.699.8255 or Scan to: BMSorders@crosstex.com

Preferred Dealer: _