

LIMITED TIME PROMOTIONAL OFFERS!

*Mix & Match: C-Files • K-Files • Flexicut Files Hedstroms • Reamers Gates Gliddens Drills Peeso Reamer Drills

Purchase 25 packs*, get

5 Packs FREE!

Promo Code RHK

Purchase 10 packs*, get

2 Packs FREE

Promo Code RHL


Purchase 5 packs*, get

1 Pack FREE!

Promo Code RHM

See back for special Aha! offers!

Continued GROWTH helps you achieve goals whether it is through education, training or practice. Roydent is committed to growing with you. That's why throughout 2018 and beyond, we are focused on delivering educational tools to help you grow your practice. From clinical tips to detailed product highlights and training, you can rely on Roydent to focus on growth.



	1140	
FILE	TIP	DESIGN
C-FILES	Non-Cutting tip	Designed for initial access and breaking through calcification.
FLEXICUT FILES	Non-Cutting tip	Triangular blade for increased flexibility. Designed for full instrumentation.
K-FILES	Cutting tip	Designed for full instrumentation.
NITI FILES	Non-Cutting tip	Conforms to the curvature of the canal during instrumentation.
Also available: Hedstroms, Reamers, Mid-Sized K-Files.		







Introducing the Aha![™] Family of Impression Materials

We understand how frustrating it can be to THINK you have the perfect impression only to have molds chip, break, or tear in processing. With Roydent's new Aha! family of impression materials, you can feel confident knowing you will get a solid, durable impression that will hold up and stay true. You get a great mold AND a happy patient, the FIRST time. Now that truly IS a great impression!





- Can be used for two-step or one-step putty-wash impression technique hydrophilic properties, high tear strength, dimensional accuracy, and a high resistance to permanent deformation.
- · High initial hydrophilicity.
- Working time: 2 min, Working time at 35°C/95°F: 1 min, Intraoral set time: 2 min, Total set time: 4 min

Aha! Regular Set Bite Registration Material NEW SKU: K11932R

- · Syringeable addition-curing elastomeric bite registration materials
- · Used for making accurate occlusal registrations.
- Shore-hardness D 43
- · Working time: 30 sec Intraoral set time: 1 min 30 sec Total set time: 2 minn

Aha! Fast Set **Bite Registration Material** NEW SKU: K11961R

- Syringeable addition-curing elastomeric bite registration materials
- Used for making accurate occlusal registrations.
- · Can be used for applications which require a hard bite registration silicone
- Shore-hardness D 43
- · Working time: 15 sec. Intraoral set time: 45 sec, Total set time: 1 min



Aha! Putty Fast Set NEW SKU: K11121R

- Standard kneadable putty
- Can be used as two-step putty wash in conjunction with Aha! VPS X-Light
- · Can be used as one-step putty wash
- · Working time: 2 min, Intraoral set time: 2 min. Total set time: 4 min

SPECIAL INTRODUCTORY OFFER: Buy 2, Get 1 FREE!

Includes all Aha! registration and VPS materials | code RW1



The boutique brand you've come to love & trust.

roydent.com | 1.800.992.7767 | 🛐 🛅 🛅





