

# **2021** Spotlight Products

JANUARY - MARCH

## Biodentine® Macrive Biosilicate Technology

#### **Bioactive Dentin Substitute**

- Resin-free formulation = no shrinkage
- Non-toxic to pulp cells, so it can be used in direct pulp caps
- Its high alkaline, pH of 12, will help to arrest any residual caries left behind

#01C0600 15 capsules & 15 single-dose pipettes • #01C0605 5 capsules & 5 single-dose pipettes

### BioRoot™ RCS MACTIVE PLOS IN ACTIVE PLOS INCLUDENT

#### **Bioceramic Root Canal Sealer**

- Pure mineral tri-calcium silicate
- · Uses cold single cone or cold lateral condensation
- · Outstanding adhesion to dentin and gutta-percha points

#01E0300 35 application pack with new double sided spoon



### NEW

## R.T.R+ (Resorbable Tissue Replacement)

#### **Biphasic Synthetic Bone Substitute**

2 formulations: 40/60 and 80/20 (substitute 80/20 for original RTR)

• Infrabony defects • Peri-implant defects • Sinus lift • Ridge augmentation • Cystic cavities

#01S0540 40% B-TCP | 60% Hydroxyapatite • #01S0530 80% B-TCP | 20% Hydroxyapatite

# **R.T.R Membrane** (Resorbable Tissue Replacement)

#### Absorbable Collagen Membrane

- Protect periodontal ligament development from fast growing soft tissues
- Maintain bone grafting material to the region to be reconstructed
- Support wound healing & guide bone regeneration

#01S0510 (15 x 20mm), 1 per box • #01S0520 (20 x 30mm), 1 per box



# **Ultra Safety Plus**

#### **Safety Needle System**

- · No recapping necessary
- Bevel indicator assists in orientating the bevel to the bone
- With transparent barrel to allow visualization of aspiration

**Buy 3 + Get 1**\* (While supplies last. Excludes Ultra Saftety handles)

**#01N4252** 25g Long (red) **#01N4271** 27g Short (orange) **#01N4272** 27g Long (yellow) **#01N4301** 30g Short (blue) **#01N3000** Safety Plus Syringe Handle **#01N3400** Pre-Sterilized Single Use Handles - 50 individually packed

# ERGONOMIC LIGHTWEIGHT CORDLESS

DENTAPENBYSEPTODONT.COM # 0 1 N 6 0 0 0





## **Septocaine®**

(Articaine Hydrochloride 4% and Epinephrine)

- Predictable duration of anesthesia: up to 60 minutes for infiltration injections and 120 minutes for nerve blocks
- Fast onset (1 to 9 minutes) to avoid delay before effective analgesia is achieved

**#01A1400** 1:100,000 Injection (Articaine Hydrochloride 4% and Epinephrine) **#01A1200** 1:200,000 Injection (Articaine Hydrochloride 4% and Epinephrine)

## **DycloPro™ Topical Anesthetic**

(Dyclonine Hydrochloride Topical Solution 0.5%)

A topical anesthetic useful to suppress the gag reflex, can also be applied to the gums prior to scaling (prophylaxis)

#99510 30g: Multiple doses in each bottle



## Septoject® Standard

Sterile single use dental needle

High grade surgical stainless steel tubing to reduce risk of breakage

#01N1252 25g long (Red) • #01N1271 27g short (Orange)
#01N1272 27g long (Yellow) • #01N1300 30g x-short (Purple) • #01N1301 30g short (Blue)

## **OraVerse**™(Phentolamine Mesylate)

The First and Only Local Anesthetic Reversal Agent
Allows you to bring patients "back to normal" sensation twice as fast

#99210 1 blister of 10 x 1.7ml glass cartridges



# DENTAPEN ELECTRONIC SYRINGE

DENTAPENBYSEPTODONT.COM



videos, CE webinars,

Download
The Septodont App
gives instant access to
product literature,